

**DE**

Dieses Medizinprodukt ist unsteril. Vor der Anwendung muss das Medizinprodukt sterilisiert werden. Die Sterilisation hat ausschließlich durch einen validierten Sterilisationsprozess mit Ethylenoxid (EO) gemäß den geltenden Bestimmungen der ISO 11135, ISO 11137-1 und ISO 10993-7 in der jeweils gültigen Version zu erfolgen. Andere Sterilisationsverfahren sind ausgeschlossen. Die Anwendung des Medizinproduktes ist nur qualifizierten Ärzten und befugtem Personal gestattet.









OP-Abdeckungen mit Klebefläche und/oder Inzisionsfolie können beim Patienten Hautirritationen hervorrufen. Die Hautempfindlichkeit ist durch den Anwender zu prüfen.






**EN**

This medical device is non-sterile. The medical device must be sterilised before use. Sterilisation must be carried out exclusively using a validated sterilisation process with ethylene oxide (EO) in accordance with the applicable provisions of ISO 11135, ISO 11137-1 and ISO 10993-7, as amended. Other sterilisation processes are excluded. Only qualified surgeons and authorised personal are permitted to use the medical device.



Surgical drapes with an adhesive surface and/or incision foil can cause skin irritation in patients. Skin sensitivity must be checked by the user.

						
DE	CE-Konformitätskennzeichen	Hersteller	Nicht steril	Herstellungsdatum	Artikelnummer	Charge
EN	CE-Marking of conformity	Manufacturer	Non-Sterile	Date of Manufacture	Article-Number	Batch Number
FR	Marque de conformité CE	Fabricant	Non stérile	Date de fabrication	Numéro d'article	Numéro de lot
IT	Marcatura di conformità CE	Produttore	Non sterile	Data di produzione	Numero articolo	Lotto

					
DE	Trocken aufbewahren	Vor Sonnenlicht schützen	Nicht wiederverwenden	Medizinprodukt	Gebrauchsanweisung beachten
EN	Keep dry	Protect from sunlight	Do not reuse	Medical Device	Consult Instructions for use
FR	Conserver au sec	Protéger de la lumière du soleil	Ne pas réutiliser	Dispositif médical	Respecter le manuel d'utilisation
IT	Conservare in luogo asciutto	Proteggere dai raggi diretti del sole	Non riutilizzare	Dispositivo medico	Osservare le istruzioni per l'uso